Wer. 1: wrzesień 2018

Nr ref. FSN: FSN-20250401-01

Nr ref. FSCA: FSCA-20250401-01

Data: 01/04/2025

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PELNAC™

Artificial Dermis: Standard Type, Fortified Type, Single Layer Type, Fenestrated Type

STERILE EO

Ethylene oxide sterilized. Do not use if package is open or damaged. Single use only. Use immediately after opening. Any portions unused after opening the package should be discarded. Do not re-sterilize. Store in a dry place. Avoid exposure to high temperatures and protect from freezing (1°C / 34°F - 30°C / 86°F). Expiry date is indicated on the outer packaging.

Description

PELNAC™ consists of two layers basically; a porcine tendon-derived atelocollagen sponge layer and silicone film. The fortified type also contains a non-adhesive silicone gauze (TREX™) mesh to reinforce the silicone film. The single layer type consists of a collagen sponge layer only. The fenestrated type has many slits into the fortified type structure to ease the drainage of exudates.

Mode of Action

When applied to full-thickness skin defects PELNAC™ encourages the formation of new dermis-like tissue by invasion of fibroblasts into the atelocollagen sponge matrix. PELNAC™ has also been shown to suppress the contraction of the wound that accompanies healing.

Intended Use

PELNAC™ is intended to be used in the following applications:

- Full-thickness burns
- Traumatic skin defect wounds
- Skin defects after tumors or nevus removal
- Site of skin flap extraction.

Indications

PELNAC™ is suitable for use in full-thickness skin defect wounds including full-thickness burns, traumatic skin defect wounds, skin defects after tumor or nevus removal and sites of skin flap extraction.

Contraindications

- PELNAC™ may exacerbate conditions in patients showing sensitivity to porcine-derived products (such as insulin), or silicone materials.
- PELNAC™ may increase infection in patients showing a sudden rise in body temperature and who appear to be showing signs of infection during the use of PELNAC™.
- Do not use in patients with a history of hypersensitivity to proteins of animal origin.
- Do not use in infected wound sites.

Precautions

- Caution should be exercised in patients susceptible to such allergic symptoms as bronchial asthma or urticaria.
- PELNAC™ has no antibacterial activity and care must be taken regarding bacterial infection. In particular, if infected wounds are present near the application site, adequate disinfection should be performed at the time of operation. If infection does occur it should be treated in accordance with local clinical practice.
- In particular, when PELNAC™ is used on a moving area such as a joint, mechanical dislodgment of PELNAC™ should be avoided.
- Discard device if mishandling has caused possible damage or contamination.
- Use PELNAC™ carefully to prevent the tear when suturing it, especially for the single layer type. Use the fortified type or the fenestrated type when the tear is expected.
- PELNAC™ should not be applied until excessive exudates, bleeding, acute swelling and infection are controlled.
- Use the fenestrated type when a lot of exudates are emitted and therefore drainage is necessary, and when a relapse of the infection in the wound in which the infection was removed is expected. [Because there is a possibility that excessive exudates separate PELNAC™ from the wound surface and that the infection relapses.]
- Use it with careful attention to prevent the bacterial intrusion, dryness and accumulation of water.

- Detach the silicone layer before the dermis-like tissue reaches the silicone layer, observing the dermis-like tissue situation from about one week after the operation. Remove the silicone layer completely surgically in case the silicone layer is involved by the dermis-like tissue, in particular when using the fenestrated type.
- Before application, thorough debridement or excision must be performed to remove any remaining necrotic tissue that may cause infection.
- If any of the following conditions occur, PELNAC™ should be removed: infection, wound colonization, sepsis, chronic inflammation (initial application of PELNAC™ may be associated with transient, mild, localized inflammation), allergic reaction, excessive redness, pain or swelling. In case of incident intra-operational or post-operational, the device should be removed and disposed based on the protocol and/or the clinical judgment of the practitioner.

Instruction for use

1. Pre-application

- (1) Confirm visually for breaches of packaging integrity of aluminum pouch prior to use.
- (2) In case the collagen sponge surface and the wound surface do not adhere closely, immerse the product in sterile physiological saline and fully saturate. When saline is infiltrated, the collagen sponge layer becomes transparent.
- (3) Perform adequate haemostasis of the wound surface; clean and disinfect it without leaving any hematoma.

2. Application

- (1) Cohesively apply one or more sheets, depending on the size and depth of the wound surface, according to the protocol of the practitioner.
- (2) Apply the collagen sponge surface to the wound surface. The silicone layer (identified by gross appearance) must be placed uppermost. Do not mistake the back and front.
- (3) Fix the product to the periphery of the wound with sutures or a surgical stapler.
- (4) When using multiple sheets, fix the sheets together in the same manner.

3. Post-application

- (1) Cover the upper surface with gauze and secure it with light pressure. Excessive pressure, dead spaces or a deviation on the wound surface prevents the invasion of capillary vessels and cells, and might result in insufficient formation of new dermis-like tissue.
- (2) Change the secondary dressing as required. Frequency of the secondary dressing change will be dependent upon the volume of exudates produced, type of dressing used and the clinician's need to inspect the wound bed for signs of infection or healing.
Note: if hematoma or excess exudates gather under the silicone film, small openings can be cut in the sheet to allow the fluid to drain.
- (3) As healing occurs, sections of the silicone film may gradually peel off and may be removed during dressing changes. Do not forcibly remove sections of the silicone film that may adhere to the wound caused by premature silicone film separation and separation of PELNAC™ from the wound bed. On inspection, if PELNAC™ is no longer covering the wound, place an additional piece of PELNAC™ over the wound.

Suturing, stapling and dressings should be applied for PELNAC™ according to the protocol and/or the clinical judgment of the practitioner.

If using chemicals (e.g. alcohol), be careful not to contact with PELNAC™ during surgery.

Disposal

Product to be disposed according to institutional procedures.

Possible Risks/Adverse Events

Possible adverse events are infection (and inflammation, fever, rash, toxic shock syndrome, infectious disease resulting from it), chronic swelling, excessive redness, excessive pain, excessive expansion, hypertrophic scar, hematoma, seroma, allergic reaction and graft failure, common complications in patients with full-thickness skin defects.

The use of more than following device area of PELNAC™ may have ethylene oxide sterilization residues (ethylene oxide and ethylene chlorohydrin) that exceed the tolerable exposure (TE) levels.

Product type	Product area not exceeding TE
Standard type	6,295 cm ²
Fortified type	7,710 cm ²
Fenestrated type	
Single layer type	27,625 cm ²
Single layer type thin thickness (1.5 mm)	65,529 cm ²

PELNAC™ may have ethylene oxide sterilization residues (ethylene oxide and ethylene chlorohydrin) that exceed the tolerable exposure (TE) levels for neonates, premature neonates and infants. Therefore, do not use PELNAC™ in excess of the area shown in the following table for neonates, premature neonates and infants.

Product Type	Product area not exceeding TE	
	Neonates weighing 3.5 kg	Premature neonates weighing 1 kg
Standard type	314 cm ²	89 cm ²
Fortified type	385 cm ²	109 cm ²
Fenestrated type		
Single layer type	1,381 cm ²	393 cm ²
Single layer type thin thickness (1.5 mm)	3,276 cm ²	933 cm ²

Reference: Area of 3L size product; 480 cm², LL; 288 cm², L; 98.4 cm², M; 73.8 cm², S; 49.2 cm², SS; 24 cm², 3S; 12 cm²

Product Availability

Standard type PELNAC™

Code	Size	Dimensions (mm)	Availability (sheet / box)
PN-R20030	4S	20 × 30 × 3	1
PN-R40030	3S	40 × 30 × 3	1
PN-R40060	SS	40 × 60 × 3	1
PN-R82060	S	82 × 60 × 3	1
PN-R82090	M	82 × 90 × 3	1
PN-R82120	L	82 × 120 × 3	1
PN-R120240	LL	120 × 240 × 3	1
PN-R200240	3L	200 × 240 × 3	1

Fortified type PELNAC™

Code	Size	Dimensions (mm)	Availability (sheet / box)
PN-F20030	4S	20 × 30 × 3	1
PN-F40030	3S	40 × 30 × 3	1
PN-F40060	SS	40 × 60 × 3	1
PN-F82060	S	82 × 60 × 3	1
PN-F82090	M	82 × 90 × 3	1
PN-F82120	L	82 × 120 × 3	1
PN-F120240	LL	120 × 240 × 3	1
PN-F200240	3L	200 × 240 × 3	1

Single layer type PELNAC™ (3mm thickness)

Code	Size	Dimensions (mm)	Availability (sheet / box)
PN-S20030	4S	20 × 30 × 3	1
PN-S40030	3S	40 × 30 × 3	1
PN-S40060	SS	40 × 60 × 3	1
PN-S82060	S	82 × 60 × 3	1
PN-S82090	M	82 × 90 × 3	1
PN-S82120	L	82 × 120 × 3	1
PN-S120240	LL	120 × 240 × 3	1
PN-S200240	3L	200 × 240 × 3	1

Single layer type PELNAC™ (1.5mm thickness)

Code	Size	Dimensions (mm)	Availability (sheet / box)
PN-S20030T	4S	20 × 30 × 1.5	1
PN-S40030T	3S	40 × 30 × 1.5	1
PN-S40060T	SS	40 × 60 × 1.5	1
PN-S82060T	S	82 × 60 × 1.5	1
PN-S82090T	M	82 × 90 × 1.5	1
PN-S82120T	L	82 × 120 × 1.5	1
PN-S120240T	LL	120 × 240 × 1.5	1
PN-S200240T	3L	200 × 240 × 1.5	1

Fenestrated type PELNAC™

Code	Size	Dimensions (mm)	Availability (sheet / box)
PN-D20030	4S	20 × 30 × 3	1
PN-D40030	3S	40 × 30 × 3	1
PN-D40060	SS	40 × 60 × 3	1
PN-D82060	S	82 × 60 × 3	1
PN-D82090	M	82 × 90 × 3	1
PN-D82120	L	82 × 120 × 3	1
PN-D120240	LL	120 × 240 × 3	1
PN-D200240	3L	200 × 240 × 3	1

Symbols

 Do not re-use	 Do not re-sterilize	 Caution
 Catalogue number	 Batch code	 Temperature limit
 Keep dry	 Keep away from sunlight	 Sterilized using ethylene oxide
 Size (length × width × thickness)	 Use-by date: yyyy-mm-dd (year/month/day)	 Do not use if package is damaged
 Manufacturer	 Authorized representative in the European Community	 Authorized representative in Switzerland
 Importer in Switzerland		

**Manufacturer:**

GUNZE LIMITED

46, Natsumegaichi, Aono, Ayabe, Kyoto 623-8513, JAPAN

Phone: +81-773-42-8035

**Authorized representative in the European Community:**

Emergo Europe B.V.

Westervoortsedijk 60, 6827 AT Arnhem, The Netherlands

Phone: +31-70-345-8570

Email: EmergoVigilance@ul.com

UK Responsible Person:

Emergo Consulting (UK) Limited

c/o Cr360 – UL International, Compass House, Vision Park Histon

Cambridge CB24 9BZ, United Kingdom

Phone: +44-1223-772-671

Email: UKRPvigilance@ul.com

**Authorized representative in Switzerland:**

MedEnvoy Switzerland

Gotthardstrasse 28, 6302 Zug, Switzerland

Phone: +41-41-562-0142

Email: EmergoVigilance@ul.com

**Importer in Switzerland:**

TapMed Swiss AG

Gumprechtstr. 33, 6376 Emmetten, Switzerland

Phone: +41-41-520-6111

Email: Info@tapmed-swiss.ch

INSTRUCTIONS FOR USE FOR:

PELNAC™

English

Deutsch

Français

Italiano

Español

Nederlands

Português

Polski

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PNEU-240425